

SENSICA


Sensifirm training

Sensifirm




What is cellulite ?

- Cellulite is the herniation of subcutaneous fat within fibrous connective tissue that manifests as skin dimpling and nodularity, mainly affecting upper thigh, buttocks and abdomen of women.
- It occurs in approximately 85-98%(*) of post-pubescent women.
- The cause for cellulite appearance can be related to: female anatomy and hormones, unwanted vascular changes resulting in a decrease in blood flow, and chronic inflammation leading to adipocytes (fat cells) lysis and skin atrophy.
- A high-stress lifestyle have also been associated with the development of cellulite.





Radio Frequency

- RF (Radio Frequency) is a type of electromagnetic energy that causes increased rotation of water molecules, leading to efficient heating of biological tissue. The thermal effect of this localized electrothermolysis process induces collagen and elastin production.
- RF energy, by penetrating deeper than laser radiation, is a superior method, not only for treatment of facial rhytides and laxity, but also for body tightening
- The RF waves reach deep into your skin, where they produce a thermal effect.
- The massage movement allows the RF waves to reach a wider range of skin layers, enabling more efficient and even tissue heating.
- Moreover, the heat induces production of new-fresh collagen, leading to accumulation of the effect over time, enabling long-lasting results.
- The depth of penetration of RF energy is inversely proportional to the frequency. Consequently, lower frequencies of RF are able to penetrate more deeply.



LIPOTHERM™ RF Technology: How It Works

Revolutionizing the traditional approach of non-invasive cellulite reduction and body contouring treatments, Sensica has developed a proprietary technology: LIPOTHERM RF™ (Radio Frequency). In a unique dual action process, while the electrodes of the device massage your skin, radio waves reach deep into the skin layers. This creates a deep thermal effect that activates collagen production and decreases the volume of adipocytes (fat cells).



RF combined with gentle vibrations enhances local stimulation, activating collagen production and decreasing the volume of fat cells. The result is visibly less cellulite, reduced body circumference and a firmer body.

[Sensifirm animation](#)





MY CELLULITE TREATMENT

A NON-INVASIVE HOME-USE TREATMENT THAT
SHRINKS FAT CELLS

- Suitable for all skin types and complexions
- Ergonomic - Shaped to imitate natural body curvature for maximum treatment comfort
- Easy to use
- Suitable for large treatment areas
- FDA cleared
- 4 electrodes
- 2 temperature safety sensors
- 3 energy levels
- Smart- even on unseen areas the vibrations will tell you that the device is operating

Sensifirm



Benefits

- Smoother and tighter-looking body
- Improves the appearance of cellulite
- Firm the skin and improve its overall look
- Shrinks fat cells
- Reduces Body Circumference(*) (avg. 2-2.5 cm)
- Stimulates collagen generation
- Stimulate cell's metabolism
- Enhanced microcirculation and lymphatic drainage thanks to vibrations assisted by a manual massage



(*) Clinical study: RD 12049 E0 Sensifirm Safety and Efficacy Study- A Five Months Report

RPL™ Before & After

Before




After



SAFETY


- LIPO THERM RF™ double-safety mechanism:
- Unique sensors constantly measure your skin's temperature and adjust RF energy delivery accordingly.
- Sensors verify that the electrodes are in full contact with the skin (at least two electrodes- rear or front). When no contact is detected RF delivery is deactivated.

Temperature sensors



Intended Use


- Sensifirm is intended for women who desire to reduce the appearance of cellulite and achieve a visibly firmer body.
- Perform a Preliminary test 48 hours prior to full first treatment!
- The desired areas is to be treated once a week for 8 consecutive weeks.
- 10–20 minutes per area (depending on size of the area).
- Protocol to be followed by maintenance treatments as needed.




How To Use

4 Quick & Easy steps to Cellulite Reduction

1. After a preliminary test - 48 hours prior to full treatment
2. Begin treatment at lowest energy level and increase only if it feels comfortable
3. Make sure that the device is in full contact with the skin during the whole treatment (green light)
4. Move continuously the device in slow circular or back and forth motions




Make sure the skin
is clean and dry



Apply Base Gel



Select RF level



Start Treatment

High User satisfaction

96%


Showed improvement in cellulite appearance after 3 months

89%

Self reported an improvement in the skin's texture, tightening and cellulite appearance

96%

Reported feeling an improvement in their skin's appearance



72%

Showed improvement in cellulite appearance already after 4 treatments

25%

Decrease in cellulite severity was measured 3 month after treatment

0%

Side effects



Sensilift



Sensimatch



Sensilight



Sensideep



Sensifirm