

# Sensilift

Your professional anti-aging treatment

PRODUCT NAME: Sensilift

PRODUCT REFERENCE: ST100EU04SNS

EAN: 5060445341279

DEEE HT :        DEEE TTC :  
ECOLOGIC



## Consumer Benefits

- Reduces wrinkles on: face, hands and décolletage
- Stimulates Collagen remodeling and Collagen production
- Stimulates Elastin production
- Improve Texture and feeling of firmness
- Visible Results after 1 treatment
- Firm and lift your skin
- Enhance microcirculation
- Enhance the appearance of skin's radiance.
- Non-invasive
- Immediate and visible results after one treatment
- Treats multiple skin layers simultaneously
- Safe, double safety mechanism

## USP vs Competitors

- Unique patented technology: Dynamic RF™
- Based on clinically proven RF technology for wrinkles reduction
- Combines facial massage for multiple layers treatments
- FDA cleared
- Two dynamic electrodes
- 2 sensors to monitor skin contact
- 2 sensors to adjust temperature for a more personalized a safe treatment
- FDA cleared
- Quick, 5-minute treatment once a week
- Ergonomic Design
- Professional technology
- Easy to use

## In the box:

- Sensilift device
- Base Gel
- Warranty card
- External power supply
- Instructions for Use & quickguide
- Cotton bag

## TECHNICAL DESCRIPTION:

140 gr.

Powered by bipolar RF at  $10 \pm 0.5$  [Watt]

**Frequency:**  $1 \pm 0.05$  [MHz]

**Mechanical massage Speed:** 50-100 RPM

**Temperature sensors:** RF, heating is automatically deactivated when skin temperature reaches  $40$  [°C]/ $104$ [F].

The heating is reactivated when temperature decreases.

**External power supply:** Input  $100-240$  [VAC],  $50-60$  [Hz], Output  $9$ [VDC],  $1.5$ [A]

Qty per master carton: 8 devices	Length cm	Width cm	Height cm	Weight Kg
Box size	22.5	22.5	10	0.99
Master dimensions	50	51	27	9