

# Sensilight MINI

Your permanent hair reduction



PRODUCT NAME : Sensilight Mini

PRODUCT REFERENCE :  
Sensilight Mini 300 - IPL030EU03SNS

EAN :  
Sensilight Mini 300 - 5060445341514

DEEE HT : DEEE TTC :  
ECOLOGIC

## CONSUMER'S BENEFIT

- Permanent hair reduction for body and face
- Professional RPL™ - Reactive Pulsed Light with skin tone and contact sensors for a safe treatment
- Fast- 1 flash every 3 seconds - full leg treatment in less than 8 minutes
- Big treatment window - 3 cm<sup>2</sup>
- Suitable for skin type 1 to 5 and most hair colours (except white grey and red)
- Up to 75% hair reduction after only 7 treatments
- Easy to use
- 2 year warranty

## USP AGAINST COMPETITORS

- Unique RPL™ technology for a safe, personalized treatment
- Compact design
- High capacity lamp for a life time of treatments
- Emitted energy: 5 J/cm<sup>2</sup>

## TECHNICAL DESCRIPTION :

Sensilight Mini device

## IN THE BOX :

Power cord  
User manual  
Warranty card

## LOGISTIC DATA

Qty per master box: 4 devices	Length cm	Width cm	Height cm	Weight Kg
Box size	18.5	18.5	10	0.69
Master dimensions	39	38	23	5.5