

# SensilightPRO

Your permanent hair reduction

## PRODUCT REFERENCE

EU IPL300EU04SNS

UK IPL300UK04SNS

## EAN

EU 5060445341439

UK 5060445341446

Unlimited  
flashes



## RPL™ TECHNOLOGY

Sensilight PRO for permanent hair reduction uses RPL™, Reactive Pulsed Light, an innovative patented technology to remove unwanted hair on the face and body. RPL™ utilizes the clinically proven IPL technology along with a combination of unique skin condition sensors for your complete safety and ease of use. IPL based technology, works by transferring light through the skin's surface, which is absorbed by the melanin in the hair shaft. The light then turns to heat energy, disabling the hair follicle and preventing regrowth.

**With two different sized treatment windows along with unlimited flashes and cordless functionality, the Sensilight Pro will give you a flexible, efficient, and effective treatment.**

## PRODUCT DETAILS

- Cordless
- 1 flash every ~1 second
- 600 flashes per charge
- Skin tone and contact sensor
- Energy recommendation button
- Precision attachment 2cm<sup>2</sup> for smaller areas
- Large treatment window 4.5cm<sup>2</sup> for large areas
- 2 operational modes: Pulse and Automatic Glide

hair  
colours



skin  
tones



## WHAT'S IN THE BOX?

- Sensilight Pro device
- External power supply
- Precision attachment 2cm<sup>2</sup>
- Instructions for Use
- Quick Start Guide
- Warranty card

## LOGISTIC DATA

Qty per master box: 4 devices	Length (cm)	Width (cm)	Height (cm)	Weight (Kg)
Box size	26	22.2	11	1.3
Master dimensions	50	30	28	5.6

SENSICA